

Visual Aid for Cartridge Inspection

Contents

1 Purpose	2
2 Scope	2
3 Reference Documents	2
4 Definitions	2
5 Responsibilities	2
6 Procedure	3
6.1. Minor Visual Defects	3
6.1.1. Cracked Damage on body of the cartridge.	3
6.1.2. Damaged or cracked lid.	4
6.1.3. Broken Lid Hinge.	5
6.1.4. Leakage from Cartridge.	5
6.1.5. Liquid observed under the film seal.	6
6.1.6. Crystallization on lid.	7
6.1.7. Crystallization on outside of cart.	8
6.1.8. Lid Overwelding.	9
6.1.9. Cracked/Broken chimneys.	10
6.1.10. Film delamination.	11
6.1.11. Poor Film Seal.	12
6.1.12. Incomplete Chimney Melt/White Crescents around the ports	13
6.1.13. Beads, Retain balls or debris in plunger hole.	14
6.1.14. Retain balls, beads, or debris in Chambers 2,3,5 and 8.	15
6.1.15. Damaged or loose Prefilter.	16
6.1.16. Upside-down prefilter	17
6.2. Major Visual Defects	18
6.2.1. Bent or Broken reaction tube	18
6.2.2. Absence of reagents in assay/product specific Chambers	18
6.2.3. Missing or Loose Large Retain Ball	19
6.2.4. Misaligned Valve Body	20
6.2.5. Chipped, Damaged or Deformed Valve Body	22
6.2.6. Missing or improper label placement	23
6.2.7. Label misprint, unreadable, faded or Unscannable Barcode	25
6.2.8. Punch in label.	25
6.2.9. Incorrect components	26
6.2.10. Missing components	27
6.3. Post-run cartridge inspection (QC testing only)	29

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 **2:54:33 PM**

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.3.1. PCR tube film fuse artifact	29
7 Document Revision History.....	31

1 Purpose

This document is for use with D8148, D16199 and D12547. It provides a visual aid for defects found on completed cartridges.

2 Scope

This document pertains to visual inspection that occurs during Manufacturing and In-Process Testing.

3 Reference Documents

D8148 Reagents On Board Assembly Line (ROBAL) Operating Instructions

D16199 ROBAL II Operating Instructions

D36739 NFC Standalone Operating Instructions

D36415 USPS Standalone Operating Instructions

D12547 ROBAL Product Testing Procedure

4 Definitions

Major Visual Defect- The defect has a likelihood to cause a major functional failure in the cartridge, Erroneous-But-Believable (EBB) result or misinformation to the end users.

Minor Visual Defect- The defect has a likelihood to cause a minor functional failure in the cartridge, Non-Determinate (ND) result or complaints regarding to the appearance of the cartridge.

5 Responsibilities

Manufacturing and In-Process Testing are responsible for performing the activities described in this procedure as well as maintaining this procedure.



Visual Aid for Cartridge Inspection

6 Procedure

6.1. Minor Visual Defects

6.1.1. Cracked Damage on body of the cartridge.

Table 1 – Cracked Damage on body of the cartridge

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 1 – Cracked Damage on the cartridge body.</p> 	<p>Acceptable Figure 1 – No Evidence of Damage on body of cartridge.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.







Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.2. Damaged or cracked lid.

Table 2 – Damaged or cracked lid

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 2. Single crack at any part of the back lid, potential for lid to leak.</p> 	<p>Acceptable Figure 2 – No evidence of a single crack on the back of the lid.</p> 
<p>Unacceptable Figure 3a - Multi-cracked lid at the back of the lid, lid particulate is breaking off from lid and multi-crack on the side of the lid, potential for lid particulate to break off and cause leaks. All multi-cracked lids at the back or sides of the lid are unacceptable</p> 	<p>Acceptable Figure 3a - No evidence of multi-cracks on the lid. Single crack on the side of the lid. No signs of leakage</p> 
<p>Unacceptable Figure 4b - Multi-cracked lid at the front of the lid, lid particulate is breaking off from lid and potential for lid particulate to break off and cause leaks.</p> 	<p>Acceptable Figure 4b - No evidence of multi-cracks on the front of the lid. No signs of leakage</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.



Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection



6.1.3. Broken Lid Hinge.

Table 3 – Broken Lid Hinge

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 5 – Broken hinge.</p> 	<p>Acceptable Figure 5 – No broken hinge.</p> 

6.1.4. Leakage from Cartridge.

Table 4 – Leakage from Cartridge

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 6 – Evidence of leakage (foot leak, lid leak).</p>  <p><i>*(Information for IPT) In cases where cartridge(s) leaked onto the rest of the cartridges in the sample bag, locate, mark, and count the leaky cartridge(s) as a Minor Visual Failure(s) and replace the remaining cartridges in the time point.</i></p>	<p>Acceptable Figure 6 – No evidence of leakage (no leak from the foot or lid).</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.



Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.5. Liquid observed under the film seal.

Table 5 – Liquid observed under the film seal

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 7 – Liquid observed under the film seal.</p> 	<p>Acceptable Figure 7 – No liquid observed under the film seal.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.





Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.6. Crystallization on lid.

Table 6 – Crystallization on lid

Unacceptable/Minor Visual	Acceptable
<p data-bbox="191 422 764 485">Unacceptable Figure 8 – Moderate crystallization observed.</p>  <p data-bbox="180 1100 802 1247">Unacceptable Figure 9 - Any presence of liquid on lid surface, any presence of crystallization or liquid on or touching the plunger hole, moderate crystallization on more than 1 chimney or crystallization on more than 5 % of lid surface.</p> 	<p data-bbox="849 422 1458 485">Acceptable Figure 8 – Light crystallization observed through the film seal.</p>  <p data-bbox="857 947 1468 1037">Acceptable Figure 9 - Thin crystallization streaks on or under the lid, light crystallization on max 1 chimney or crystallization on max 5 % of lid surface.</p> 

Visual Aid for Cartridge Inspection

6.1.7. Crystallization on outside of cart.

Table 7 – Crystallization on outside of cart

Unacceptable/Minor Visual	Acceptable
<p data-bbox="180 422 773 478">Unacceptable Figure 10 – Moderate crystallization observed.</p> 	<p data-bbox="846 422 1455 451">Acceptable Figure 10 – No crystallization observed.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.





Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.8. Lid Overwelding.

Table 8 – Lid Overwelding, e.g., plastic strings

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 11 – Lid Overwelding with Plastic Hair.</p> 	<p>Acceptable Figure 11 – Stress whitening, small lid lip melting</p> 
<p>Unacceptable Figure 12 - Lid Overwelding, Plastic Strings, Damage to Upper and Front Lid</p> 	<p>Acceptable Figure 12 - Stress whitening, small lid lip melting.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.






Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.9. Cracked/Broken chimneys.

Table 9 – Cracked/Broken chimneys.

Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 13 – Cracked/Broken chimneys.</p> 	<p>Acceptable Figure 13 – Normal chimney appearance.</p> 
<p>Unacceptable Figure 14 - Missing Chimney</p> 	
<p>Unacceptable Figure 15 - Chimney moved sideways / mashed chimney.</p> 	
<p>Unacceptable Figure 16 - Chimney moved under the lid</p> 	

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.



Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.1.10. Film delamination.

Table 10 - Film delamination.

Unacceptable/Minor Visual	Acceptable
<p data-bbox="168 426 784 485">Unacceptable Figure 17 - Film delamination (flap) on any chimney or lid.</p>  <p data-bbox="168 898 802 978">NOTE: Whether delamination is acceptable or not depends on the affected assay and the functional impact, see D12547 for more information</p>	<p data-bbox="873 426 1430 485">Acceptable Figure 14 - No delamination on any chimney or lid.</p> 

Visual Aid for Cartridge Inspection

6.1.11. Poor Film Seal

NOTE: To differentiate lift/void from glare, view top-down then tilt 10–20°. Glare moves with angle; a true lift/void tends to remain.

NOTE: Poor outer seal (outer film seal) is referred to as POFS in PQC context.

Table 11A – Poor inner seal.



Unacceptable/Minor Visual	Acceptable
<p>Unacceptable Figure 18 - Poor inner seal</p> 	<p>Acceptable Figure 15 –Good inner and outer seal</p> 

Table 11B – Poor outer seal (Outer film seal / POFS).

Unacceptable/Minor Visual	Acceptable
---------------------------	------------

Visual Aid for Cartridge Inspection

Unacceptable Figure 19 - Poor outer seal



Acceptable Figure 16 –Good inner and outer seal



6.1.12. Incomplete Chimney Melt/White Crescents around the ports

Table 12 – Incomplete Chimney Melt

Unacceptable/Minor Visual	Acceptable
---------------------------	------------

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

Unacceptable Figure 19 – Incomplete Chimney Melt



Acceptable Figure 15 – Complete Chimney melt



Unacceptable Figure 20 – White Crescents around the ports.

6.1.13. Beads, Retain balls or debris in plunger hole.

Table 13 – Beads, Retain balls or debris in plunger hole

Unacceptable/Minor Visual	Acceptable
---------------------------	------------

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

Unacceptable Figure 21 – Beads or debris in plunger hole.



Acceptable Figure 16 – No beads or debris in plunger hole.



6.1.14. Retain balls, beads, or debris in Chambers 2,3,5 and 8.

Table 14 – Retain balls, beads, or debris in Chambers 2,3,5 and 8.

Unacceptable/Minor Visual

Unacceptable Figure 22 – Retain balls, beads, or debris in sample chamber (Chamber 3) or liquid reagent chambers (Chamber 2, 5, and 8).



Acceptable



Acceptable Figure 17 – Reagent chambers and sample chambers (Chambers 2, 3, 5, and 8) are clear of beads and debris.



Visual Aid for Cartridge Inspection

6.1.15. Damaged or loose Prefilter.





Table 15 – Damaged or loose Prefilter

Unacceptable/Minor Visual	Acceptable
<p data-bbox="159 447 795 506">Unacceptable Figure 23 – Damaged Prefilter (e.g., out of place, cracked, loose).</p> 	<p data-bbox="837 447 1474 506">Acceptable Figure 18 – No Damage (e.g., Intact Prefilter).</p> 

Visual Aid for Cartridge Inspection

6.1.16. Upside-down prefilter

Table 16 - Upside-down prefilter

Unacceptable/Minor Visual	Acceptable
<p data-bbox="193 424 760 483">Unacceptable Figure 204 – Wrong orientation in chamber</p>  <p data-bbox="170 1037 784 1096">Unacceptable Figure 25 – Standalone filter in wrong orientation.</p> 	<p data-bbox="834 424 1466 453">Acceptable Figure 19 - Correct orientation in chamber</p>  <p data-bbox="857 1068 1451 1127">Acceptable Figure 20 – Standalone filter in correct orientation.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM



Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.2. Major Visual Defects



6.2.1. Bent or Broken reaction tube

Table 17 - Bent or Broken reaction tube

Unacceptable/Major Visual	Acceptable
<p>Unacceptable Figure 26 - Damage (e.g., bent, broken reaction tube or broken reaction tube tip).</p> 	<p>Acceptable Figure 21 - No damage (e.g., reaction tube intact and reaction tube tip intact).</p> 

6.2.2. Absence of reagents in assay/product specific Chambers

Table 18 - Absence of reagents in assay/product specific Chambers

Unacceptable/Major Visual	Acceptable
<p>Unacceptable Figure 27 - Absence of liquids on board reagents, if applicable (e.g., Chamber 2, 5 and 8).</p> 	<p>Acceptable Figure 22 - Presence of liquids on board reagents, if applicable (e.g., Chamber 2, 5 and 8).</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.



Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.2.3. Missing or Loose Large Retain Ball

Table 19 – Missing or Loose Large Retain Ball (USPS Only)

Unacceptable/Major Visual	Acceptable
<p data-bbox="180 422 776 483">Unacceptable Figure 28 – Loose Large Retain Ball visible at film seal</p> 	<p data-bbox="844 422 1456 483">Acceptable Figure 23 – Large Retain Ball is properly inserted</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.





Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.2.4. Misaligned Valve Body

Table 20 - Misaligned Valve Body (major visual)

Unacceptable/Major Visual	Acceptable
<p data-bbox="168 489 779 552">Unacceptable Figure 29 – Valve body is misaligned. (This is a more easily determined failing example).</p>  <p data-bbox="191 1129 779 1192">Unacceptable Figure 29a – Arrow pointing outside the H.</p> 	<p data-bbox="873 489 1433 552">Acceptable Figure 24a - Valve body is correctly aligned. Arrow pointing at the H.</p>  <p data-bbox="883 1123 1450 1186">Acceptable Figure 24b – Valve body is correctly aligned. Arrow pointing, but borderline, at the H.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

Unacceptable Figure 29b – Arrow pointing outside the H.



Acceptable Figure 24c – Valve body is correctly aligned. Arrow pointing, but borderline, at the H.



Acceptable Figure 24d – Valve body is correctly aligned. Arrow pointing, but borderline, at the H.





***Note:** The arrow needs to be pointing within the H. The arrow pointing at the line, or the line edge of the H is acceptable.

Visual Aid for Cartridge Inspection

6.2.5. Chipped, Damaged or Deformed Valve Body

Table 21 - Chipped, Damaged or Deformed Valve Body

Unacceptable/Major Visual	Acceptable
<p data-bbox="245 422 711 485">Unacceptable Figure 30 - Valve body is chipped/damaged/deformed.</p> 	<p data-bbox="837 422 1468 453">Acceptable Figure 25- Valve body shows no damage.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.



Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.2.6. Missing or improper label placement

Table 22 - Missing or improper label placement (non-NFC cartridge)

Unacceptable/Major Visual (non-NFC cartridge)	Acceptable (non-NFC cartridge)
<p>Unacceptable Figure 31 – Unacceptable cartridge label placement (e.g., missing label, improper placement).</p> 	<p>Acceptable Figure 26 - Acceptable cartridge label (e.g., label present, proper placement).</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.



Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection



Table 23 - Missing or improper label placement (NFC cartridge)

Unacceptable/Major Visual (NFC cartridge)	Acceptable (NFC cartridge)
<p data-bbox="131 390 776 449">Unacceptable Figure 32 - Touching blue lid/foot and/or hanging off to the side.</p>  <p data-bbox="139 1619 768 1675">Not a failure at IPT (see section 5.6.5.2 in D12547 for more information)</p>	<p data-bbox="878 405 1482 434">Acceptable Figure 27 - Label present and centered.</p> 

Visual Aid for Cartridge Inspection



6.2.7. Label misprint, unreadable, faded or Unscannable Barcode

Table 24 - Label misprint, unreadable, faded or Unscannable Barcode

Unacceptable/Major Visual	Acceptable
<p>Unacceptable Figure 33 - Cartridge info incorrect, faded, or barcode cannot be scanned (e.g., illegible, unscannable).</p> 	<p>Acceptable Figure 28 - Barcode is scannable and cartridge information (e.g., Batch, Lot, Cartridge S/N) is clear, legible, and displays the right information.</p> 

6.2.8. Punch in label

Table 25 - Punch in label (For NFC products ONLY)

Unacceptable/Major Visual (NFC cartridge)	Acceptable (NFC cartridge)
<p>Unacceptable Figure 34 - Punch in label.</p>  <p>Classified as a barcode failure at IPT (see section 5.6.5.2 in D12547 for more information)</p>	<p>Acceptable Figure 29 - No Punch.</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary


Visual Aid for Cartridge Inspection

6.2.9. Incorrect components

Table 26 - Incorrect components (wrong lid, PCR-tube, or VBA-tube).

If a cartridge with incorrect components is discovered, it is unacceptable/major visual

Lids



PCR-tubes



VBA-tubes



User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM


Confidential - Company Proprietary

Visual Aid for Cartridge Inspection


6.2.10. Missing components

Table 27 - Missing components (missing plunger, missing pre-filter, missing funnel, etc.).

If a cartridge with missing components is discovered, it is unacceptable/major visual.





The VBA assembly includes the plunger and a cartridge without a plunger will not run.



Visual Aid for Cartridge Inspection

6.2.11 Missing Center Perforation

Table 28 – Missing Center Perforation on Top Film Seal (USPS Only)

Unacceptable/Major Visual	Acceptable
<p data-bbox="167 457 787 514">Unacceptable Figure 35 – Missing Center perforation on top film seal</p> 	<p data-bbox="841 457 1461 514">Acceptable Figure 30 – Center perforation on top film seal</p> <p data-bbox="870 516 1432 548">Note: Image is magnified for better image quality</p> 

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary

Visual Aid for Cartridge Inspection

6.3. Post-run cartridge inspection (QC testing only)

6.3.1. PCR tube film fuse artifact

6.3.1.1. Also sometimes referred to as PCR tube doughnut/donut.

6.3.1.2. This inspection applies ONLY to post-run inspection of QC functionally tested cartridges if potentially qualifying for a re-test per D12547 due to PCR tube film fuse artifact.

6.3.1.2.1. This inspection applies only to cartridges that underwent seal test AND functional testing and resulted in specific errors, see details in D12547.

6.3.1.2.2. If the post-run cartridge is not available for inspection or inspection result is inconclusive, the replicate may NOT be re-tested. A cartridge which may have been physically damaged post-run is considered not available for inspection.

6.3.1.2.3. PCR tube film fuse is not considered a visual failure but only used to motivate meeting criteria for re-test per D12547.

6.3.1.3. A cartridge is considered to have a PCR tube film fuse if meeting the following criteria:

6.3.1.3.1. The PCR tube is filled with liquid. The fused area must be surrounded by liquid.

6.3.1.3.2. The PCR tube film for a part of the PCR tube window is fused/stuck together with no liquid between the two film areas. This fused area must be at least 1mm in size (1 dimension) and may have different shapes (commonly diamond or oval shaped similar to the shape of the PCR tube window) and is usually located in the middle of the PCR tube. Refer to Table 29 for examples.

6.3.1.3.3. Consult with FI if guidance is needed.





6.3.1.4. The following observations do not qualify for PCR tube film fuse: Air bubble in PCR tube, precipitates in PCR liquid, wavy PCR tube film. Refer to Table 29 for examples.

6.3.1.4.1. Importantly, a PCR tube film fuse will not move while air bubbles or precipitates could move upon e.g., turning or gently tapping the cartridge, or gently touching of the PCR tube (which should be generally avoided and if, should be done under utmost care not to damage the PCR tube). Air bubbles may be observed additionally to PCR tube film fuses.

6.3.1.4.2. White/milky shadows may sometimes be observed on PCR films post-run. These are not considered PCR tube film fuses if no actual fuse is observed and there is still liquid between the two film pieces.

Visual Aid for Cartridge Inspection

Table 29 – Post-run PCR tube film fuse inspection

NOT qualifying for PCR tube film fuse	Qualifying for PCR tube film fuse
<p>Figure 36 – Visual findings that do NOT indicate a PCR tube film fuses:</p> <ul style="list-style-type: none"> PCR tube not liquid filled (observed best by comparing the light refraction of the tube to a non-run PCR tube empty cartridge and a post-run passing cartridge) Air bubble (without PCR tube fuse)  <ul style="list-style-type: none"> PCR liquid precipitates  <ul style="list-style-type: none"> Wavy PCR films 	<p>Figure 31 – Examples of PCR tube film fuses:</p>  <ul style="list-style-type: none"> Fused area must be at least 1mm in size (1 dimension) PCR tube film fuse post-run may be combined with air bubbles (often multiple smaller ones)

Visual Aid for Cartridge Inspection

7 Document Revision History

Rev	Effective	Description of Changes
Y.1	Current	Correcting accidental text formatting in Rev. Y (red and underlined text in section 6.3)
Y	9/24/2024	Addition of Section 6.3. Post-run cartridge inspection (QC testing only): PCR tube film fuse artifact to help with implementation of re-test for PQC failures caused by PCR tube film fuse in D12547.
W	3/21/2024	Addition of Sec 6.1.12 with the Table 12 including Figure: 15,19, and 20 to help associates to easier identification of acceptable vs unacceptable incomplete chimney melt/white crescents around the ports.
V	03/08/2024	Added Table 27 – Missing Center Perforation on Top Film Seal *USPS Only
U	11/30/2023	Replaced Figure 28 with Figure 28a. Replaced Figure 23 with Figure 23a. Added Figures 28b, 23b, 23c and 23d to 6.2.4 to help associates to easier determine acceptable vs. unacceptable VBA misalignment.

User of a printed copy of this document is responsible for using the current revision and discarding expired or obsolete revisions.

Current released version of this document resides in Agile Workplace.

Cepheid Confidential Date Printed: 1/12/2026 2:54:33 PM

Confidential - Company Proprietary